

IHE Integration Statement		Date:	February 2008
Vendor	Product Name	Version	
Carestream Health Inc.	Kodak DryView 6800 Laser Imager	1.02	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server	None	
Mammography Image (MAMMO)	Print Server	None	
Internet address for Carestream Health IHE information: http://www.carestreamhealth.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestreamhealth.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document