

|                                                                                                                                                            |                                                                             |                                           |          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|-------------------------------------------|----------|
| IHE Integration Statement                                                                                                                                  |                                                                             | Date:                                     | Sep 2019 |
| Vendor                                                                                                                                                     | Product Name                                                                | Version                                   |          |
| Carestream Health                                                                                                                                          | Carestream DRX Revolution Mobile System and Carestream DRX Evolution System | 1.0, Build 1.5                            |          |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |                                                                             |                                           |          |
| Integration Profiles implemented                                                                                                                           | Actors Implemented                                                          | Options Implemented                       |          |
| Scheduled Workflow (SWF)                                                                                                                                   | Acquisition Modality                                                        | None                                      |          |
| Patient Information Reconciliation (PIR)                                                                                                                   | Acquisition Modality                                                        | None                                      |          |
| Consistent Presentation of Images (CPI)                                                                                                                    | Acquisition Modality                                                        | None                                      |          |
| Radiation Exposure Monitoring (REM)                                                                                                                        | Acquisition Modality                                                        | None                                      |          |
| Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html                                                                 |                                                                             |                                           |          |
| Links to Standards Conformance Statements for the Implementation                                                                                           |                                                                             |                                           |          |
| DICOM                                                                                                                                                      | http://www.carestream.com/dicom.html                                        |                                           |          |
| Links to general information on IHE                                                                                                                        |                                                                             |                                           |          |
| North America:<br>http://www.ihe.net                                                                                                                       | Europe:<br>http://www.ihe-europe.net                                        | Japan:<br>http://www.jira-net.or.jp/ihe-j |          |

End of Document

---