

| IHE Integration Statement                                                                                                                                  |                                                                                                                                                                                                                                                                                                                                                          | Date:                                                                                   | Dec 2019 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|----------|
| Vendor                                                                                                                                                     | Product Name                                                                                                                                                                                                                                                                                                                                             | Version                                                                                 |          |
| Carestream Health                                                                                                                                          | CARESTREAM DRX-Ascend System<br>Quantum Q-Rad-DIGITAL DRX-Series                                                                                                                                                                                                                                                                                         | 1.0, Build 1.5                                                                          |          |
| Carestream Health                                                                                                                                          | CARESTREAM Compass OTC system                                                                                                                                                                                                                                                                                                                            | 1, Build 1.6                                                                            |          |
| Carestream Health                                                                                                                                          | CARESTREAM DRX-Ascend System<br>Quantum Q-Rad-DIGITAL DRX-Series<br>CARESTREAM DRX-1 System<br>CARESTREAM Mobile X-ray System<br>CARESTREAM DRX-Revolution Mobile X-ray System<br>CARESTREAM DRX-Revolution Nano Mobile X-ray System<br>CARESTREAM Motion Mobile X-ray System<br>CARESTREAM DRX-Evolution System<br>CARESTREAM DRX-Evolution Plus System | 1, Build 1.7                                                                            |          |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |                                                                                                                                                                                                                                                                                                                                                          |                                                                                         |          |
| Integration Profiles implemented                                                                                                                           | Actors Implemented                                                                                                                                                                                                                                                                                                                                       | Options Implemented                                                                     |          |
| Scheduled Workflow (SWF)                                                                                                                                   | Acquisition Modality                                                                                                                                                                                                                                                                                                                                     | None                                                                                    |          |
| Patient Information Reconciliation (PIR)                                                                                                                   | Acquisition Modality                                                                                                                                                                                                                                                                                                                                     | None                                                                                    |          |
| Consistent Presentation of Images (CPI)                                                                                                                    | Acquisition Modality                                                                                                                                                                                                                                                                                                                                     | None                                                                                    |          |
| Radiation Exposure Monitoring (REM)                                                                                                                        | Acquisition Modality                                                                                                                                                                                                                                                                                                                                     | None                                                                                    |          |
| Internet address for Carestream Health IHE information: <a href="http://www.carestream.com/ihe.html">http://www.carestream.com/ihe.html</a>                |                                                                                                                                                                                                                                                                                                                                                          |                                                                                         |          |
| Links to Standards Conformance Statements for the Implementation                                                                                           |                                                                                                                                                                                                                                                                                                                                                          |                                                                                         |          |
| DICOM                                                                                                                                                      | <a href="http://www.carestream.com/dicom.html">http://www.carestream.com/dicom.html</a>                                                                                                                                                                                                                                                                  |                                                                                         |          |
| Links to general information on IHE                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                          |                                                                                         |          |
| North America:<br><a href="http://www.ihe.net">http://www.ihe.net</a>                                                                                      | Europe:<br><a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>                                                                                                                                                                                                                                                                             | Japan:<br><a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a> |          |

End of Document