

IHE Integration Statement		Date:	November 2007
Vendor	Product Name	Version	
Carestream Health	Kodak POC CR Systems	3.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Composer	None	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
Internet address for Kodak IHE information: http://www.carestreamhealth.com/go/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestreamhealth.com/go/dicom		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

Pub No.