

IHE Integration Statement		Date:	January 2007
Vendor	Product Name	Version	
Carestream Health	Kodak DirectView DR Systems	4.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Internet address for Carestream Health IHE information: http://www.carestreamhealth.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestreamhealth.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document