

IHE Integration Statement		Date:	15 Aug 2003
Vendor	Product Name	Version	
Eastman Kodak Company	Kodak DryView 8200 Laser Imager with Kodak PACS Link Systems	5.2	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server	None	
Basic Security (SEC)	Secure Node	None	
Internet address for Kodak IHE information: www.kodak.com/go/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	www.kodak.com/go/dicom		
Links to general information on IHE			
North America: www.rsna.org/IHE	Europe: www.ihe-europe.org	Japan: www.jira-net.or.jp/ihe-j	