


# Dental Film




Behind the box you see today  
is a long heritage of expertise

The same pioneering  
**quality** and **reliability**,  
only the name has changed

Early 20<sup>th</sup> century



Today



### Kodak, the original dental film for more than a century

A true original, Kodak film changed the way the world of dentistry works. Our experience has allowed us to continuously push the standards of image quality further.

We were once known as Kodak Dental Systems, but now we are a part of Carestream Dental. Nothing about the film itself has changed. We will continue our role as a dynamic pioneer contributing to the further development of the world of dental film and digital imaging.

### The Carestream Dental advantage

- Comfortable and patient friendly
- High, consistent batch to batch quality
- A broad range of products available everywhere
- Comprehensive dealer network
- Easy to open packaging

# Kodak's heritage is now part of Carestream Dental's unique approach

## Carestream Dental delivers the difference to your practice

Carestream Dental is dedicated to enhancing dentists' everyday workflow with the aid of solutions designed to perfectly reflect their needs. We define this unique approach to the oral health care market as the Carestream Dental Factor. Exceptional film and digital image quality is a cornerstone of our dedication; it is the tool that enables dentists to make the right diagnosis and determine the optimal treatment.



### [Workflow integration]

Our systems and software are specifically designed to enhance care and to fit seamlessly into your busy practice. This allows you to spend less time on administrative tasks and more time providing patient care.

### [Humanized technology]



We keep our technology simple for you and your patients. Clear communication is the key to obtaining patient confidence and optimizing patient care.

### [Diagnostic excellence]


To make faster, more accurate diagnoses, detail is everything. Our film, diagnostic systems, digital imaging and software do more than just help you identify a problem; they show you where to take action.

The **best quality...**

...in a **broad range**



**INTRAORAL FILM**



**INSIGHT**

The standard for intraoral film. With INSIGHT dental film, you'll get a true F-speed film that requires up to 60% less radiation per exposure.



**Ultra-speed**

A D-speed film that generates the high-quality, detailed images you need to make accurate diagnoses.



**E-Speed**

Blending quality and speed, E-speed provides brilliant image clarity with as little as one-half the exposure time required for traditional D-speed intraoral films.



**D-Speed**

A cost-effective, high-quality solution if you need the technique settings of a traditional dental radiographic system. It delivers fine detail, high contrast and outstanding processing and handling characteristics.

**EXTRAORAL FILM**



**T-MAT G/RA**

This dental film provides high-contrast, detailed images of intervening tooth structures while retaining good soft tissue visibility when used with LANEX regular or medium screens.



Would you like to **know more?**  
Visit **[carestream.com](https://www.carestream.com)**